

Lien vers project Proposal :

<https://docs.google.com/document/d/1-zwgUUmswFLc1M6NGqPrPLSwb89mPJKqguFEqY7JhRM/edit>

## Lab Notebook :

25/01/2017

Choice of a project : how *Caenorhabditis elegans* react to an electric field.


- Does the electric field change the way *C. elegans* move ?
- Is *C. elegans* able to orientate through the field ?

Remarks from Tamara :

- We should work with electrophoresis and turn the gel where *C. elegans* are
- Print in plastic, above the gel, channels and use them to count the worms per direction.
- We can't use a voltmeter it's not precise and adapted for the project, we should use an electric field sensor made with antenna and transistor.


Choice of the set up:

- Average speed = 400  $\mu\text{m/s}$
- Size of *C. elegans* = around 1 mm
- Straight microchannels (5 cm long and 80  $\mu\text{m}$  deep) with varying widths of 2 mm, 1 mm, 500  $\mu\text{m}$ , 300  $\mu\text{m}$ , and 150  $\mu\text{m}$  with electrodes instrumented in their reservoirs (extracted from *Electrotaxis of Caenorhabditis elegans in a microfluidic environment*)
- About the microscope :



### Observations / Quantifications :

- Speed (measuring the distance travelled / time) => ImageJ
- Type of movements (?)
- Orientation (counting how many *C. elegans* are in a channel, cf. draw 2)



- About the sensor : quantify the intensity, the orientation, the current
- About the electric field : AC/DC, Voltage, Current, Signal shape (square, triangle, ...)

### Pre-experiment :


- We checked the singer microscope to do our experiments.
- We will use a 96 wells plate box (without the wells of course) to realise our experiments, it is possible to see 8cm long, yet we only need to see 5cm long.
- We can see the *C. elegans* quite well and can record with the video. We can move the video to track the *C. elegans*.
- We found by calculation that since a plate radius is 4,5cm and we generally must pour 20mL of agar to fill it, then the depth of the medium is 3mm.
- Optimal frequency : 1mHz to 100mHz ( between 100mHz and 1kHz it's not working very efficiently)
- Shape of signal is not supposed to have influence on the *C. elegans* mvmt
- 20°C
- Circle = 2,5cm away from each electrode

### Media for the “electrophoresis”:

“The agar gel was composed of: de-ionized water, 2% of Bacto-Agar, glycerol (3.7 mL of glycerol 60% for 1 L), NaCl (0.250 mmol/L) as previously described in [16]. The gel was cast by pouring a first layer of agar and adding a PDMS (polydimethylsiloxane) block onto it so that it will shape the future cavity where nematodes will crawl (6×8 cm). A second layer of gel was then poured around the PDMS block. Once solidified the PDMS block was removed. The resulting agar pad was then placed in an electrophoresis box filled with a buffer. It was composed of de-ionized water, glycerol (3.7 mL of glycerol 60% for 1 L) and NaCl (0.250 mmol/L) as previously described in [16]”



To prepare the plates where *C. elegans* will be observed, we prepare media for 10 plates. We measure that 60 mL of water is needed to fill the plates with 3 mm height which . We take into account that the structure printed will also reduce the volume by half.


### About the set up :



First design of the channels, in the center we put the worms and then let them go in one of the eight channels. The purpose is to see both speed and direction of the worms depending the frequency/amplitude/orientation of the electrical field.

\*There are more structure to come\*






This chip will be used to see if the average speed of *C. elegans* changes with frequency/amplitude of the electrical field.

**Settings : frequency - amplitude - AC/DC**

**Measurements : speed by measuring distance and time**

**What is fixed : orientation of the field**




This chip will be used to see if the orientation of lines of the electrical field induces the worms to change their direction.

**Settings : orientation of the field**

**Measurements : number of individuals in each channel**

**What is fixed : frequency and amplitude**





This chip will be used to see the sensitivity of the *C. elegans* to the location of the electric field source.


**Settings : orientation of the field**

**What is fixed : frequency - amplitude - AC/DC**

**Measurements : speed by measuring distance and time**




These are the **two controls** used : the first, the negative one, prevent worms from moving and the second one makes them go straight .




The chip is put inside the square plate and recovered partially by the medium.

Remarks : In each experiment, the current must be fixed around 50 mA.



We expect the lines of field to be as shown in this picture (*from wikipedia*)



## Source of power :

Digital Dual channel DDS Function Signal Generator MHS-2300A 5MHZ



from aliexpress

Here is the datasheet :

<http://blog.pennybuying.com/download/fb/FB170A.pdf>

## Medium preparation

Before the experiment we prepare the medium for the “electrophoresis”. We use a variant of the recipe see before.

For 500ml of media:

osmosis water : 500ml

Glycerol 60% : 1,85ml

NaCl : 0,125mmol → 0,025g for 500ml

Agar (2%) → 10g

## Pre-experiment:

We poured this medium, on our square plate, with our models in wood with the channel of 2mm and 1mm.

After a battery of tests with the models, we concluded that the medium can go in all the channels, therefore we decided to start the experiment with the C.elegans on Monday.

# 29/01/17

We decided to use a new electronic sensor : a Movuino, able to detect a magnetic field and therefore an electric field.

The python code is almost finished.

# 30/01/17

## Experiment:

During our experiment we try to observe the *C.elegans* with the microscope for yeast (micromanipulator), but this microscope is too precise, because we can't see the totality of the channel. We choose to change the experiment. We just cut a little piece of agar and we put this piece under the binocular magnifying glass, who can take the totality of the piece, and see the *C.elegans*.

We change the power device too, at the beginning we use Digital Dual channel DDS Function Signal Generator MHS-2300A 5MHZ, but he doesn't work, so we try to use a battery of 9V, but it doesn't work too, and at the end we choose use the device [...].

But during this day we could manipulate the *C.elegans*, begin to observe their behaviour, movement etc and adapt the way to observe them.

# 31/01/2017

The totality of our expectation about the reaction of *C.elegans* in function of electric field are dead.

We try different things with the agars, the power supply and the microscope and nothing work.

We try to take one by one the *C.elegans* and transfer a part of 1,5 cm X 2cm (the maximum size of viewing of our microscope), but we don't have the material to do this right, and we can't manipulate them. We try to 'glue' the two part of agar, the virgin part of agar without *C.elegans* and the other with the 'colony' of *C.elegans* to permit to *C.elegans* to go in the virgin part. And after 'glue' put the electric field with the negative point in the virgin part, wait 5minute, and see how many *C.elegans* are in the ex-virgin part. But we can't really 'glue' the two part and the *C.elegans* can't passed.

We try to put few *C.elegans* in the virgin plate to scrapping the first plate thanks to a trombone.

This technique work, but at the end we can't capture the picture, and we don't observe any special behaviour.

To have data, we try to change the experiment, and do the same things, but with *Pyrocystis fusiformis* a bioluminescent algae.

[...]

At the end of the day, we don't have any data or new experiment, and more we don't have new experiment.

Therefore we decided to test the electric field on all the organisms that we could use in the lab.

- One of the protocol involved dinoflagella, organisms that react to movement by luminescence. We found literature about how electromagnetic fields would influence the luminescence of this organisms.  
Therefore we tried to put a plate containing 10mL of media and organisms, and put the plate on a vortex.  
The idea was to induce movement with the vortex and apply an electric field to influence the luminescence. Sadly, despite this experiment being amazing, the luminescence could hardly be observed, nor measured.
- We also tried to open and close carnivorous plants, *dionaea*, by inducing an electric field between the two parts of it, with the idea that since the signal is an electric signal, we could force the signal into the plant.  
Unlucky us it didn't work quite as planned, however it still might be due to a low intensity of the power supply.
- We then tried to do the same experiment of the *C. elegans*, with *Chironomus plumosus* and *Artemia salina*. In order to do this we placed them in water, each in different plates of course. We then applied the electric field and took videos of each plates.

Our results are based on the study of the videos of *Chironomus plumosus*, compared to the study of the sensitivity of *movuino*.